ID N° 20032 - Contact Angle Measurement - M4^s

Executed on Feb 17, 2023, by M. Breitwieser

Wasser - Standard-Wilhelmy-Platte (5mm-Scan)

• Wilhelmy-Platte (DIN) / Wasser : (20,0°C, 1,85', 5,0mm, 0,0787mm/s - static)

$\theta_{C,s} = 3,9 \pm 0,92^\circ$ CAH 5° , ${}^aE_s = 72,524 \pm 0,041$ mN/m ${}^aH_s = -0,367$ mN/m, $B_{\%,s} = 99,7\%$

Report

1. ===== Introduction =====

Principle: The M4 contact angle measurement using the *modified Wilhelmy method* is based on force measurements on a cuboid or cylindrical solid body brought into contact with a liquid surface. The contact line between the solid, the liquid surface, and the gas phase is called the three-phase boundary line (triple line), and the angle at which the liquid surface meets the three-phase boundary line is the contact angle θ (theta). The relative position of the three-phase boundary line-above, on, or below the surrounding liquid level-reflects the surface energy characteristics and is indicated by the resulting curvature of the meniscus. The weight of this meniscus curvature is determined gravimetrically, taking buoyancy forces into account, and is related to the length of the three-phase boundary line. With perfect wetting ($\theta=0^\circ$), this specific force reaches a maximum value, which corresponds to the surface tension of the liquid; with poor wettability ($\theta>90^\circ$), however, this wetting force becomes negative. Energetic interactions between the solid and liquid phases are thus made visible through the contact angle. When the solid body is immersed, the *triple line scans the solid surface*, yielding the advancing contact angle (θ_A), and when it is withdrawn, the receding contact angle (θ_R) is obtained - measured by force of wetting. The difference ' $\theta_A-\theta_R$ ' is called the contact angle hysteresis (CAH).

Physically: The Wilhelmy equation $\cos\theta_{(A|R)} = {}^aF / (p \cdot \gamma)$ describes the relationship between the contact angle θ , the perimeter of the solid p , the wetting force aF , and the surface tension of the liquid γ . Adhesion tension aE is defined according ${}^aE_{A,R} = m_{meniscus} \cdot g / p = {}^aF / p = \gamma \cdot \cos\theta = \gamma_s - \gamma_{sl}$

Practically: In this measurement, after the position of contact of liquid and solid has been determined, the 'Wilhelmy-Platte (DIN)' test specimen is moved into and out of the Wasser liquid by incremental movements. This allows the determination of *static contact angles*. Furthermore, a defined force stabilization can be maintained before each subsequent movement. This enables *static contact angle* measurement based on the patent DE 4412405, which provides the often required equilibrium measurements in accordance with physical principles. This control technology yields a stabilization time for levelling and for each of the selected 34 measuring positions, the respective contact angle $\theta_{A,s}$ or $\theta_{R,s}$ is calculated independently.

Information: The IMETER M4 method description, available at '<https://imeter.de/?view=article&id=46>', provides information on procedures and calculations.

2. ===== Static Contact Angle =====

Wilhelmy-Platte (DIN), cubic plate 20x0,2mm

5,0mm Wilhelmy-Platte (DIN) \ Wasser, $\gamma = 72,74$ mN/m $\vartheta = 20,0^\circ$ C

Contact Angle, CA θ_C	$3,9 \pm 0,92^\circ$	$\theta_{M,s} = 4^\circ$	
Contact Angle Hysteresis, CAH _s	5°		
	Advancing $\theta_{A,s}$	Receding $\theta_{R,s}$	Transition - mm
Contact Angles, θ_s	$6,0^\circ$	$1,5^\circ$	$0^\circ_{adv.} \Rightarrow 0^\circ_{rec.}$
Rated measurements $n_{mm-range}$	20 0,180 - 3,991mm	14 2,787 - 0,180mm	2 4,993 - 4,918mm
avg. Platform speed \bar{v}_z [mm/s]	$0,0787 \pm 0,0230$	$-0,0651 \pm 0,372$	
avg. Step distance Δh_0 [mm]	$0,200 \pm 0$	$-0,125 \pm 0,0612$	
avg. equilibration Time $\bar{t}_{eq.}$ [s]	$1,94 \pm 0,277$	$1,60 \pm 0,241$	

Static Contact Angle Measurement: The equilibrium Contact Angle θ_C is calculated from advancing and receding values according to the formula of Tadmor/Chibowski; the average value of advancing and receding Contact Angle is given by θ_M . The avg. equilibration Time $t_{eq.}$ is the mean of the time between positioning and acquiring the wetting force. The criterion of force equilibrium is defined in the utilized IMPro. Diagramm(s) below presenting these results and further details in a graphical way.

- Diagramm 1.2: 'Contact Angle²¹' -- Data, Calculation and Results


Diagram 'Contact Angle²¹' gives a summary on this contact angle measurement: The small triangular marks show the individual contact angle values depending on the immersion depth of the sample. At '0' on the x-axis, the flatface front of the sample touches the surface of the liquid, to the right the depth of immersion increases. So 'immersion' indicates the nominal position of the three-phase boundary ('Triple Line') on the sample surface. Red triangles show the contact angles when the fluid advances over the sample surface, Green triangles when retracting (receding). The direction reversal shows a small step, that is the 'contact angle hysteresis' (CAH). The two contact angles and the extent of hysteresis are determined from the graph by plotting linear regressions over the ranges. The advancing angle θ_A is determined by the regression equation for the immersion depth '0'; the receding angle θ_R is determined by extrapolation to the maximum

immersion depth (thus any errors due to adherence of remaining liquid are eliminated). The standard deviation of the regression curves is indicated by the width of the hatching, which also graphically indicates the data range of the regression. Further marks indicate the determined angles (red, green) and the hysteresis range (yellow) above the Y-axis, the main result, the equilibrium contact angle θ_c is shown in blue.

Adhesion Force of the equilibrated Triple Line

Energy of Adhesion, ${}^a\bar{E}_s = ({}^aE_A + {}^aE_R) / 2$ **72,524 ± 0,041 mN/m**

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
3. ===== Details on the Measurement and Setup =====

Sample: 'Wilhelmy-Platte (DIN)', **Form:** cubic plate, face dimension 20 x 0,2 mm.

Test Liquid: 'Wasser' at 19,988°C: density 0,998208 g/mL, surface tension 72,74 mN/m, viscosity 1,002 mPa·s; Capillary length 2,73 mm.

Suspension/Fixation: by operator; measured weight at start 1,2118 g.

Vessel: tempered double wall measuring cell, height 140 mm, inclosing sample and liquid at coherent temperature and atmosphere; Surface area 1452 mm².

Conditioning: during 0,9 min in the measuring cell above the surface of the liquid.


Sample Weight: at start 1,5305 g, after conditioning 1,5308 g, after the measurement 1,5337 g.

Gain in weight during conditioning 0,0003 g (+0,20‰) and after measurement 0,0032 g (+0,21%).

Measurement algorithm: static CA-measurement, acquisition of equilibrated values, stepwise movement. Maximum immersion of 5,000 mm and force equilibration at the inflection point for 2,2 s. Duration for immersion 1,1 min, for emersion 45 s.

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
'Temperature & Events²' -- Survey on Temperature and Time

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- Diagramm 4.3: 'RawData²ⁿ' -- Acquired weights during immersion and emersion (raw data and immersion corrected lengths)

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4. ===== Table of raw data and results =====

The table below provides the data for each contact angle measuring position in this run ($n=79$). - Within a row the column t lists the time of the CA-measurement. $t_{eq.}$ is the duration reach the equilibrated force after positioning (static measurement). The nominal distance between the sample bottom flatface and the liquid surface is given by h_0 , whilst $h_{corr.}$ shows the corrected immersion depth. With W_{RAW} the stabilized total weight of the sample, sample holder, wetting force and buoyancy is written. The $W_{corr.}$ -values* are former weighing values when the weight of the sample and its holder is subtracted. After buoyancy- correction, aF is the force

exerted on the triple line. The Energy of Adhesion ${}^aE_{A,R} = {}^aF/p$ is the force of adhesion per meter of the circumference (p) of the sample. θ is the Contact Angle. The last column indicates by **chr** the advancing resp. receding state; **cb** tags values of a 'capillary bridge' above the surface-level.

(Resolving static force of wetting from ${}^aF = (W_{RAW} - W_0 - W_{Buoy.}) \cdot g + F_{Buoy.air}$) * $W_{corr.} = W_{RAW} - W_0$ (with $W_0 = 2,7426$ g).

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Execution & Audit-Trail

Data created during execution of the IMPro "Wilhelmy-Sequenz, 16.02.23, 45713", type 4/4. Time Period of the Accomplishment: Feb., 17. 2023 between 16:32:23 and 16:39:38, elapsed time: seven minutes. IMPro finished as projected. The complete Report first was presented on Feb., 17.23 at 16:40. Original data was changed as reported to Audit-Log:

Start of this IMPro 16:32:46, - Backup of individual IMPro configuration, as follows:

Time Variables: _Equi_Time 0,350s, ConditioningTime 45s, StirringTime 0s.

Abs. Positions/Way: Immersion_Depth 5mm, Platform_StartPosition 87mm.

Rel. Positions/Way: FeedRateAdv 0,2mm, FeedRateRec -0,2mm, FeedRate_for_adv_rec_turn -0,075mm.

Others: _Equi_Crit 0,0002g, _Equi_loops 10x, _Equi_loopsRepeater 1x, Stirrer-Speed 0rps.

y/n: (Auto)Repeat (Auto)Repeat, Record_Air_Density -1, Temperature-Outside-Warning -1, Automatic_Thermostat_available -1.

**** CHANGE(ES) IN DATA SHEET - by user M. Breitwieser **** Day/Time: 31.05.2026 19:10:17 ****

[1. Title : from] 'das Wasser nochmal ' [to] 'Wasser - Standard-Wilhelmy-Platte (5mm-Scan)'

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The Number **20032** refers to the Recordset in the Database '*imeterData37*' where all Information can be retrieved at any time.

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